

Reliable Tools for Toxicology & Pain Management Labs

Deep-well Plates



| Color | Well volume | Well shape | Well bottom | Qty./Pack | Cat. No. |
|-------------------------------------|-------------|------------|-------------|-----------|----------|
| Deep-well plates | | | | | |
| 96-well, polypropylene | 0.3 mL | round | U-bottom | 100 | 701330 |
| 96-well, polypropylene | 0.5 mL | round | U-bottom | 48 | 701346 |
| 96-well, polypropylene | 1.1 mL | round | U-bottom | 24 | 701350 |
| 96-well, polystyrene | 1.1 mL | round | U-bottom | 32 | 701352 |
| 96-well, polypropylene, Low Profile | 1.2 mL | round | U-bottom | 50 | 701340 |
| 96-well, polypropylene | 2.2 mL | square | U-bottom | 24 | 701354 |
| 96-well, polypropylene | 2.2 mL | square | V-bottom | 50 | DW2200 |
| 96-well, polypropylene, Sterile | 2.2 mL | square | V-bottom | 50 | DW2250 |
| 384-well, polypropylene | 0.3 mL | square | V-bottom | 48 | 701355 |

Secure Plate Sealing

Mat Covers

| Description | Autoclavable (121°C for 15 min.) | Qty./Pack | Cat. No. |
|----------------------------------------------|-------------------------------------|-----------|----------|
| Pierceable Sealing Mat for 96-well | | 50 | M2200 |
| Pierceable Sealing Mat for 96-well, Sterile | | 50 | M2250 |
| Silicone cover for 384 well, deep well plate | | 50 | 701357 |
| Mat for 0.5 mL plates, PP | Yes | 50 | 701358 |
| Mat for 1.1 mL plates, mod. PE | Yes | 24 | 701360 |
| Mat for 1.2 mL plates, TPE | Yes | 50 | 701368 |
| Mat for 2.2 mL plates, EVA | No | 24 | 701362 |



Sealing Films*

| Application | Material | Qty./Pack | Cat. No. |
|-------------------------------------------|-----------------------|---------------|----------|
| Real-Time PCR and Protein Crystallization | Polyolefin, non-tacky | 100 sheets | 781391ES |
| Pierceable and Resealable | Vinyl | 100 sheets | 701396ES |
| Sealing film, aluminum | Aluminum | 1 roll of 100 | 781381 |
| Cold Storage and Automation | Aluminum | 100 sheets | 701395ES |
| DMSO Sample Storage | Polypropylene | 100 sheets | 701367 |
| Sealing Paddle for Sealing Films | | 2 paddles | 701381 |



*More available, visit www.brandtech.com

Precise, Slim, Light, & Dependable Repeating Pipette

HandyStep® S Repeating Pipette



The HandyStep **S** provides unprecedented versatility in a mechanical repeating pipette. It offers 59 different volume settings with accuracy, precision and dependability.

- Comfortable, slim design is lightweight, well-balanced and offers easy one-handed volume selection
- Flexible; accepts third party tips – no more being locked into a sole source
- Tip ejection system allows tip removal without user handling

HandyStep® touch and HandyStep® touch **S** models also available

| Description | Cat. No. | |
|--------------------------------------------------------------|-------------------------|-----------------------|
| HandyStep S Repeating Pipette | 705110 | |
| HandyStep touch Repeating Pipette | 705200 | |
| PD-Tip™ II size | Non-sterile tips | BIO-CERT Pure* |
| 0.1 mL (100 tips) | 705700 | 705730 |
| 0.5 mL (100 tips) | 705702 | 705732 |
| 1.0 mL (100 tips) | 705704 | 705734 |
| 1.25 mL (100 tips) | 705706 | 705736 |
| 2.5 mL (100 tips) | 705708 | 705738 |
| 5 mL (100 tips) | 705710 | 705740 |
| 10.0 mL (100 tips) | 705712 | 705742 |
| 12.5 mL (100 tips) | 705714 | 705744 |
| 25 mL (50 tips) | 705716 | — |
| 25 mL (25 tips) | — | 705746 |
| 50 mL (25 tips) | 705718 | 705748 |
| Multi-pack (20 tips each 0.5, 1, 1.25, 2.5, 5, 10 & 12.5 mL) | 705720 | — |

*BIO-CERT is certified sterile, and DNA-, RNase-, ATP- and endotoxin-free.

Lightweight and Rugged Pipettes

Transferpette S Adjustable Volume Single Channel Pipettes



| Capacity, μ L | Cat. No. |
|----------------------|----------|
| 0.1-1 | 705868 |
| 0.1-2.5 | 705869 |
| 0.5-10 | 705870 |
| 2-20 | 705872 |
| 5-50 | 705873 |
| 10-100 | 705874 |
| 20-200 | 705878 |
| 100-1000 | 705880 |
| 500-5000 (0.5-5 mL) | 705882 |
| 1000-10000 (1-10 mL) | 705884 |

- Simple, thumb-tip volume adjustment, even with gloves!
- Shortest stroke length among leading pipettes reduces strain
- Slim design reaches bottom of 15 mL conicals
- Comfortably fits most hands
- Accepts tips from most quality manufacturers
- Fixed volume single channels also available

Transferpette S -8/-12 Multichannel Pipettes



| Capacity, μ L | 8-Channel | 12-Channel |
|-------------------|-----------|------------|
| 0.5-10 | 705900 | 705920 |
| 5-50 | 705906 | 705926 |
| 10-100 | 705908 | 705928 |
| 20-200 | 705910 | 705930 |
| 30-300 | 705912 | 705932 |

Transferpette® S -8/-12 pipettes supplied with performance certificate, operating manual, shelf mount, two tip boxes, reagent reservoir, pack of replacement O-rings. One-year warranty.




Test Tube Racks



Polypropylene Test Tube Racks

Sturdy polypropylene racks with alphanumeric positions for easy sample identification. Racks are stackable, non-conductive and will not float in water baths. They are suitable for use in temperatures from -20°C to 90°C. Racks are supplied in packs of five. They are supplied in three pieces for easy and permanent assembly.



Test tube racks, PP

| Maximum Tube Diameter, mm | Format | Length mm | Width mm | Height mm | White Cat. No. | Blue Cat. No. | Red Cat. No. |
|---------------------------|--------|-----------|----------|-----------|----------------|---------------|--------------|
| 13 | 6 x 14 | 265 | 126 | 75 | 4340000 | 4340001 | — |
| 16 | 5 x 11 | 265 | 126 | 75 | 4340060 | 4340061 | 4340062 |
| 18 | 5 x 11 | 265 | 126 | 75 | 4340010 | 4340011 | 4340012 |
| 20 | 4 x 10 | 265 | 126 | 75 | 4340020 | 4340021 | 4340022 |
| 25 | 4 x 8 | 265 | 126 | 88 | 4340030 | 4340031 | 4340032 |
| 30 | 3 x 7 | 265 | 126 | 88 | 4340040 | 4340041 | 4340042 |

Certified Volumetric Flasks and Graduated Cylinders

BRANDTECH offers BLAUBRAND® USP volumetric instruments that comply with the Class A error limits required by the United States Pharmacopeia.

BLAUBRAND Volumetric Flask, Class A, USP, Certified, Glass Stopper

| Volume (mL) | Error limit (+/-mL) | Stopper NS | Qty./Pack | Cat. No. |
|-------------|---------------------|------------|-----------|----------|
| 5 | 0.02 | 10/19 | 2 | 36968 |
| 10 | 0.02 | 10/19 | 2 | 36973 |
| 20 | 0.02 | 10/19 | 2 | 36975 |
| 25 | 0.03 | 10/19 | 2 | 36977 |
| 50 | 0.05 | 12/21 | 2 | 36978 |
| 100 | 0.08 | 14/23 | 2 | 36979 |
| 200 | 0.1 | 14/23 | 2 | 36980 |
| 250 | 0.12 | 14/23 | 2 | 36981 |
| 500 | 0.2 | 19/26 | 2 | 36982 |
| 1000 | 0.3 | 24/29 | 1 | 36983 |
| 2000 | 0.5 | 29/32 | 1 | 36984 |

BLAUBRAND Graduated Cylinder, Tall Form, Class A, USP, Certified

| Volume (mL) | Subdivision (mL) | Error limit (+/-mL) | Height (mm) | Qty./Pack | Cat. No. |
|-------------|------------------|---------------------|-------------|-----------|----------|
| 5 | 0.1 | 0.05 | 115 | 2 | 32805 |
| 10 | 0.2 | 0.1 | 140 | 2 | 32808 |
| 25 | 0.5 | 0.17 | 170 | 2 | 32820 |
| 50 | 1 | 0.25 | 200 | 2 | 32828 |
| 100 | 1 | 0.5 | 260 | 2 | 32838 |
| 250 | 2 | 1 | 335 | 2 | 32848 |
| 500 | 5 | 2.0 | 365 | 2 | 32854 |
| 1000 | 10 | 3.0 | 465 | 1 | 32862 |
| 2000 | 20 | 6.0 | 505 | 1 | 32864 |

For more BLAUBRAND® USP glassware products see brandtech.com

