


Transferpette®S
single-channel
Nominal volume

Transferpette®S
multi-channel
Nominal volume

Pipette tips Filter tips Volume range	Transferpette®S single-channel														Transferpette®S multi-channel				
	1 µl	2.5 µl	10 µl	20 µl	25 µl	50 µl	100 µl	200 µl	500 µl	1000 µl	2 ml	5 ml	10 ml	10 µl	50 µl	100 µl	200 µl	300 µl	
0.1 – 20 µl	✓	✓ ₂₎	✓ ₂₎											✓					
0.5 – 20 µl		✓	✓											✓					
1 – 50 µl		✓ ₂₎	✓											✓					
2 – 200 µl ¹⁾				✓	✓	✓	✓	✓							✓	✓	✓	✓	
5 – 300 µl				✓	✓	✓	✓	✓							✓	✓	✓	✓	
50 – 1000 µl ¹⁾									✓	✓									
0.5 – 5 ml											✓	✓							
1 – 10 ml													✓						
0.1 – 1 µl	✓	✓ ₂₎	✓ ₂₎											✓					
0.5 – 10 µl		✓	✓											✓					
1 – 20 µl			✓											✓					
2 – 20 µl				✓	✓	✓	✓	✓							✓	✓	✓		
5 – 100 µl				✓	✓	✓	✓	✓							✓	✓	✓	✓	
5 – 200 µl								✓							✓	✓	✓	✓	
50 – 1000 µl									✓	✓									



✓ = Tip volume less than pipette's nominal volume
¹⁾ Bulk tips are yellow or blue colored; racked tips are colorless in a yellow-colored or blue-colored tip tray
²⁾ Calibration and possible adjustment needed